
# Sample BodyLogic<sup>™</sup> Reports Available on the Horizon<sup>™</sup> DXA System




## Advanced Body Composition® Report



#### Legend

- **1.** Visual image of precise location of bone, lean mass, and fat mass
- 2. Plot of % Body Fat compared to age group
- **3.** Precise Fat Mass, Lean Mass, and % Body Fat measurements of each area:
  - Left Arm
  - Right Arm
  - Torso (Trunk)
  - Abdomen (Android)
  - Hips (Gynoid)
  - Left Leg
  - Right Leg
- **4.** Calculation of Total Mass, Fat Mass, and Lean Mass to give overall % Body Fat
- 5. Estimated amount of visceral fat (the type of fat around internal organs associated with medical disorders such as metabolic syndrome, cardiovascular disease, and type 2 diabetes)

## Track % Body Fat, Fat Mass and Lean Mass values over time



### Legend

- 1. Graph of % Body Fat values over time
- 2. % Body Fat values over time
- 3. Fat Mass values over time
- 4. Lean Mass values over time
- **5.** Visual comparison over time of changes in bone, lean mass, and fat mass

MISC-04256 Rev. 001 © Hologic 2016. All rights reserved. Printed in USA. Specifications are subject to change without prior notice. Hologic, Advanced Body Composition, Horizon, The Science of Sure and associated logos are trademarks and registered trademarks of Hologic, Inc. and/or its subsidiaries, in the United States and/or other countries. All other trademarks are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to womenshealth@hologic.com.